

COSMETICS INDUSTRY

Mark 3 Moisture Analyzer

Application Examples

Fast, accurate moisture measurement with the Mark 3 allows adjustments to be made in the processing cycle in minutes rather than hours.

- Over 10,000 different materials already correlated to reference methods
- No pre-programming required on the part of the operator
- One button operation to perform a test
- Free loaners for the life of the instrument
- Two year complete parts and labor warranty
- 300 materials can be stored and recalled by name
- Hard copy results with standard internal printer
- Free applications development for the life of the instrument



TYPE MATERIAL	OVEN REFERENCE	Mark 3 RESULTS	TEST TIMES	STD DEV
Nail Polish	31.28% S	31.31 % S	5.4 Min.	.073
Nitro-base Solution	27.70% S	27.68 % S	7.0 Min	.120
Polish Remover	.273% S	0.268% S	15.3 Min.	.034
Talc	.070% M	0.069% M	2.4 Min.	.004
Styling Spritz	9.43% S	9.44 % S	5.5 Min.	.070
Anti-Bacterial Soap	43.52% M	43.497% M	5.6 Min.	.121
Body Lotion	17.51% M	17.45 % M	12.0 Min.	.160
Collagen Slurry	93.54% M	93.555% M	11.6 Min.	.030
Body Wash	21.57% S	21.62% S	11.5 Min.	.151

Table demonstrates correlation results between convection or vacuum oven and the Mark 3 moisture analyzer.

5 Orville Dr. Suite 200
Bohemia, NY 11716
www.sartorius-omnimark.com



sartorius
mechatronics

Tel (303) 876-5396
(800) 835-3211
Fax (631) 253-5263